START SAVING MONEY TODAY!

PALLADIUM SILVER ALLOY

As robust as any gold alloy.

COMPOSITION

70.9% Pd

10.0% Ag



- Only 10% Silver
- Clinically Proven
- Works With All Regular-Fusing Porcelains
- Exceptional strength for any PFM application including implants



With Argelite 71, you get a high strength alloy that can be used for both pressing and layering porcelains. The low silver content and light oxide give beautiful aesthetics, while the low density of 10.6 means more units per ounce.

Start saving money today with Argelite 71.



ARGEN'

Alloy Makers to the World

Customer Service (800) 255-5524 Technical Service (800) 255-5095

CANADA Argen Canada (566) 722-7436 Qüebec Province (514) 984-1808

www.argen.com



